

Product Datasheet

Rat ADAMTS1 ELISA Kit (Ready to Use) (orb554555)

Description	Rat A Disintegrin And Metalloproteinase With Thrombospondin 1 (ADAMTS1) Sandwich ELISA KitReady to U...
Reactivity	Rat
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q68EJ2
Sensitivity	30.5pg/mL
Entrez	79252
Expiration Date	Please enquire.