

Product Datasheet

Rat ADAMTS4 ELISA Kit (Ready to Use) (orb554556)

Description	Rat A Disintegrin And Metalloproteinase With Thrombospondin 4 (ADAMTS4) Sandwich ELISA KitReady to U...
Reactivity	Rat
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9ESP7
Sensitivity	1.15ng/mL
Entrez	66015
Expiration Date	Please enquire.