

Product Datasheet

Rat ADAMTS7 ELISA Kit (Ready to Use) (orb554558)

| | |
|--------------------------|---|
| Description | Rat A Disintegrin And Metalloproteinase With Thrombospondin 7 (ADAMTS7) Sandwich ELISA KitReady to U... |
| Reactivity | Rat |
| Range | 0.156-10ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids. |
| Sample Types | tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | Q1EHB3 |
| Sensitivity | 0.058ng/mL |
| Entrez | 315879 |
| Expiration Date | Please enquire. |