

Product Datasheet

Rat ACPP ELISA Kit (Ready to Use) (orb554564)

Description	Rat Acid Phosphatase, Prostatic (ACPP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.069ng/mL
Expiration Date	Please enquire.