

Product Datasheet

Rat ACTHR ELISA Kit (Ready to Use) (orb554581)

Description	Rat Adrenocorticotrophic Hormone Receptor (ACTHR) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	G3V848
Sensitivity	0.234ng/mL
Entrez	282839
Expiration Date	Please enquire.