

Product Datasheet

Rat AFM ELISA Kit (Ready to Use) (orb554583)

Description	Rat Afamin (AFM) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and s...
Reactivity	Rat
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P36953
Sensitivity	0.93ng/mL
Entrez	282708
Expiration Date	Please enquire.