

Product Datasheet

Rat ALDOC ELISA Kit (Ready to Use) (orb554595)

Description	Rat Aldolase C, Fructose Bisphosphate (ALDOC) Sandwich ELISA Kit Ready to Use: Pre-diluted detection ...
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P09117
Sensitivity	0.283ng/mL
Entrez	24191
Expiration Date	Please enquire.