

Product Datasheet

Rat Ab1-42 ELISA Kit (Ready to Use) (orb554618)

Description	Rat Amyloid Beta Peptide 1-42 (Ab1-42) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	6.1pg/mL
Expiration Date	Please enquire.