

Product Datasheet

Rat APOC2 ELISA Kit (Ready to Use) (orb554655)

Description	Rat Apolipoprotein C2 (APOC2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Rat
Range	39-2500pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	G3V8D4
Sensitivity	14pg/mL
Entrez	292697
Expiration Date	Please enquire.