

Product Datasheet

Rat bEP ELISA Kit (Ready to Use) (orb554696)

Description	Rat Beta-Endorphin (bEP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernatants and other biological fluids.
Sensitivity	4.35pg/mL
Expiration Date	Please enquire.