

Product Datasheet

Rat BNP ELISA Kit (Ready to Use) (orb554715)

Description	Rat Brain Natriuretic Peptide (BNP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P13205
Sensitivity	8.44pg/mL
Entrez	25105
Expiration Date	Please enquire.