

Product Datasheet

Rat C-Peptide ELISA Kit (Ready to Use) (orb554721)

Description	Rat C-Peptide Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simp...
Reactivity	Rat
Range	125-8000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	33.27pg/mL
Expiration Date	Please enquire.