

Product Datasheet

Rat CES1 ELISA Kit (Ready to Use) (orb554737)

Description	Rat Carboxylesterase 1 (CES1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q63108
Sensitivity	0.107ng/mL
Entrez	29225
Expiration Date	Please enquire.