

Product Datasheet

Rat CT1 ELISA Kit (Ready to Use) (orb554742)

Description	Rat Cardiotrophin 1 (CT1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q63086
Sensitivity	5.3pg/mL
Entrez	29201
Expiration Date	Please enquire.