

Product Datasheet

Rat CASP14 ELISA Kit (Ready to Use) (orb554750)

Description	Rat Caspase 14 (CASP14) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D3ZZ65
Sensitivity	0.103ng/mL
Entrez	299587
Expiration Date	Please enquire.