

Product Datasheet

Rat CASP7 ELISA Kit (Ready to Use) (orb554752)

Description	Rat Caspase 7 (CASP7) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O88550
Sensitivity	0.117ng/mL
Entrez	64026
Expiration Date	Please enquire.