

## Product Datasheet

### Rat CCK8 ELISA Kit (Ready to Use) (orb554774)

<b>Description</b>	Rat Cholecystokinin 8, Octapeptide (CCK8) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
<b>Reactivity</b>	Rat
<b>Range</b>	15.6-1000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
<b>Sensitivity</b>	4.45pg/mL
<b>Expiration Date</b>	Please enquire.