

Product Datasheet

Rat CHGA ELISA Kit (Ready to Use) (orb554780)

Description	Rat Chromogranin A (CHGA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q5PPG5
Sensitivity	12.3pg/mL
Entrez	24258
Expiration Date	Please enquire.