

Product Datasheet

Rat C1INH ELISA Kit (Ready to Use) (orb554814)

Description	Rat Complement 1 Inhibitor (C1INH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents fo...
Reactivity	Rat
Range	7.81-500ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6P734
Sensitivity	2.26ng/mL
Entrez	295703
Expiration Date	Please enquire.