

Product Datasheet

Rat CRHR1 ELISA Kit (Ready to Use) (orb554831)

Description	Rat Corticotropin Releasing Hormone Receptor 1 (CRHR1) Sandwich ELISA Kit Ready to Use: Pre-diluted d...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35353
Sensitivity	0.123ng/mL
Entrez	58959
Expiration Date	Please enquire.