

Product Datasheet

Rat CRHR2 ELISA Kit (Ready to Use) (orb554832)

Description	Rat Corticotropin Releasing Hormone Receptor 2 (CRHR2) Sandwich ELISA Kit Ready to Use: Pre-diluted d...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P47866
Sensitivity	0.134ng/mL
Entrez	64680
Expiration Date	Please enquire.