

Product Datasheet

Rat CTXI ELISA Kit (Ready to Use) (orb554835)

Description	Rat Cross Linked C-Telopeptide Of Type I Collagen (CTXI) Sandwich ELISA KitReady to Use: Pre-diluted...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.042ng/mL
Expiration Date	Please enquire.