

Product Datasheet

Rat CDK5 ELISA Kit (Ready to Use) (orb554841)

Description	Rat Cyclin Dependent Kinase 5 (CDK5) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q03114
Sensitivity	0.059ng/mL
Entrez	140908
Expiration Date	Please enquire.