

Product Datasheet

Rat DRD1 ELISA Kit (Ready to Use) (orb554878)

| | |
|--------------------------|---|
| Description | Rat Dopamine Receptor D1 (DRD1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for c... |
| Reactivity | Rat |
| Range | 0.156-10ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids. |
| Sample Types | tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | P18901 |
| Sensitivity | 0.062ng/mL |
| Entrez | 24316 |
| Expiration Date | Please enquire. |