

Product Datasheet

Rat DAT ELISA Kit (Ready to Use) (orb554881)

Description	Rat Dopamine Transporter (DAT) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for co...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P23977
Sensitivity	0.142ng/mL
Entrez	24898
Expiration Date	Please enquire.