

## Product Datasheet

# Rat Neutrophil Elastase ELISA Kit (Ready to Use) (orb554887)

<b>Description</b>	Rat Elastase 2, Neutrophil (Neutrophil Elastase) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
<b>Reactivity</b>	Rat
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<b>D4A488</b>
<b>Sensitivity</b>	0.069ng/mL
<b>Entrez</b>	<b>299606</b>
<b>Expiration Date</b>	Please enquire.