

Product Datasheet

Rat EPHB3 ELISA Kit (Ready to Use) (orb554903)

Description	Rat Ephrin Type B Receptor 3 (EPHB3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D3ZH39
Sensitivity	0.067ng/mL
Entrez	287989
Expiration Date	Please enquire.