

Product Datasheet

Rat FATP1 ELISA Kit (Ready to Use) (orb554932)

Description	Rat Fatty Acid Transport Protein 1 (FATP1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P97849
Sensitivity	0.071ng/mL
Entrez	94172
Expiration Date	Please enquire.