

Product Datasheet

Rat FGF10 ELISA Kit (Ready to Use) (orb554946)

Description	Rat Fibroblast Growth Factor 10 (FGF10) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Rat
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P70492
Sensitivity	33pg/mL
Entrez	25443
Expiration Date	Please enquire.