

Product Datasheet

Rat FGF6 ELISA Kit (Ready to Use) (orb554955)

Description	Rat Fibroblast Growth Factor 6 (FGF6) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8R5L5
Sensitivity	6.6pg/mL
Entrez	170700
Expiration Date	Please enquire.