

Product Datasheet

Rat gSAP ELISA Kit (Ready to Use) (orb554976)

Description	Rat Gamma-Secretase Activating Protein (gSAP) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D3ZPQ9
Sensitivity	0.113ng/mL
Entrez	311984
Expiration Date	Please enquire.