


## Product Datasheet

# Insulin Receptor Beta (phospho-Tyr1185) antibody (orb5550)

<b>Description</b>	Insulin Receptor Beta (phospho-Tyr1185) antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-Fr
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor Beta around the phosphorylation site of Tyr1185 (DI(p-Y)ET)
<b>Target</b>	Phospho-Insulin Receptor Beta (Tyr1185)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	68/152 kDa
<b>CAS Number</b>	3
<b>Uniprot ID</b>	<b>P06213</b>
<b>Expiration Date</b>	12 months from date of receipt.



Immunohistochemical analysis of formalin...