

Product Datasheet

Rat GRIA1 ELISA Kit (Ready to Use) (orb555000)

Description	Rat Glutamate Receptor, Ionotropic, AMPA 1 (GRIA1) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P19490
Sensitivity	6.1pg/mL
Entrez	50592
Expiration Date	Please enquire.