

Product Datasheet

Rat GLS ELISA Kit (Ready to Use) (orb555007)

Description	Rat Glutaminase (GLS) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P13264
Sensitivity	0.107ng/mL
Entrez	24398
Expiration Date	Please enquire.