

Product Datasheet

Rat GPX1 ELISA Kit (Ready to Use) (orb555009)

Description	Rat Glutathione Peroxidase 1 (GPX1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents f...
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P04041
Sensitivity	0.63ng/mL
Entrez	24404
Expiration Date	Please enquire.