

Product Datasheet

Rat GPC3 ELISA Kit (Ready to Use) (orb555025)

Description	Rat Glypican 3 (GPC3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P13265
Sensitivity	0.067ng/mL
Entrez	25236
Expiration Date	Please enquire.