

Product Datasheet

Rat GH ELISA Kit (Ready to Use) (orb555039)

Description	Rat Growth Hormone (GH) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01244
Sensitivity	0.054ng/mL
Entrez	24391
Expiration Date	Please enquire.