

Product Datasheet

Rat ChT ELISA Kit (Ready to Use) (orb555071)

Description	Rat High Affinity Choline Transporter (ChT) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Rat
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9GZV3
Sensitivity	28pg/mL
Entrez	60482
Expiration Date	Please enquire.