

Product Datasheet

Rat Serotonin Transporter ELISA Kit (Ready to Use) (orb555083)

Description	Rat Huntingtin (Serotonin Transporter) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	G3V9P7
Sensitivity	0.055ng/mL
Entrez	29424
Expiration Date	Please enquire.