

Product Datasheet

Rat INS ELISA Kit (Ready to Use) (orb555102)

Description	Rat Insulin (INS) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and ...
Reactivity	Rat
Range	140.62-9000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	37.3pg/mL
Expiration Date	Please enquire.