

Product Datasheet

Rat IDE ELISA Kit (Ready to Use) (orb555103)

Description	Rat Insulin Degrading Enzyme (IDE) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents fo...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35559
Sensitivity	0.129ng/mL
Entrez	25700
Expiration Date	Please enquire.