

Product Datasheet

Rat INSL3 ELISA Kit (Ready to Use) (orb555111)

Description	Rat Insulin Like Protein 3 (INSL3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Rat
Range	0.234-15ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9WUK0
Sensitivity	0.062ng/mL
Entrez	114215
Expiration Date	Please enquire.