

Product Datasheet

Rat IL17B ELISA Kit (Ready to Use) (orb555138)

Description	Rat Interleukin 17B (IL17B) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	G3V8K9
Sensitivity	6.4pg/mL
Entrez	116472
Expiration Date	Please enquire.