

Product Datasheet

Rat IAPP ELISA Kit (Ready to Use) (orb555162)

Description	Rat Islet Amyloid Polypeptide (IAPP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Rat
Range	7.81-500pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P12969
Sensitivity	2.19pg/mL
Entrez	24476
Expiration Date	Please enquire.