

Product Datasheet

Rat KISS1 ELISA Kit (Ready to Use) (orb555181)

Description	Rat Kisspeptin 1 (KISS1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenie...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q7TSB7
Sensitivity	8.47pg/mL
Entrez	289023
Expiration Date	Please enquire.