

Product Datasheet

Rat LTF ELISA Kit (Ready to Use) (orb555187)

Description	Rat Lactoferrin (LTF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D3ZAB1
Sensitivity	12.7pg/mL
Entrez	301034
Expiration Date	Please enquire.