

Product Datasheet

Rat LRP8 ELISA Kit (Ready to Use) (orb555211)

Description	Rat Low Density Lipoprotein Receptor Related Protein 8 (LRP8) Sandwich ELISA Kit Ready to Use: Pre-di...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D3ZE75
Sensitivity	13.3pg/mL
Entrez	362558
Expiration Date	Please enquire.