

Product Datasheet

Rat MIP3a ELISA Kit (Ready to Use) (orb555218)

Description	Rat Macrophage Inflammatory Protein 3 Alpha (MIP3a) Sandwich ELISA Kit Ready to Use: Pre-diluted dete...
Reactivity	Rat
Range	6.25-400pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P97884
Sensitivity	2.45pg/mL
Entrez	29538
Expiration Date	Please enquire.