

Product Datasheet

Rat MMP3 ELISA Kit (Ready to Use) (orb555234)

Description	Rat Matrix Metalloproteinase 3 (MMP3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents...
Reactivity	Rat
Range	1.25-80ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P03957
Sensitivity	0.48ng/mL
Entrez	171045
Expiration Date	Please enquire.