

Product Datasheet

Rat MTL ELISA Kit (Ready to Use) (orb555257)

Description	Rat Motilin (MTL) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and ...
Reactivity	Rat
Range	93.7-6000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A8IRIO
Sensitivity	42.5pg/mL
Entrez	126231
Expiration Date	Please enquire.